

MS3830

Wheelchair Scale



The high-capacity MS3830 wheelchair scale is portable and easily transported on rolling castors. Optional handrails can be attached for added stability when used by patients who have difficulty standing.

By deducting the weight of the wheelchair using the scale's tare function, a patient's exact weight can be accurately measured.



Castor wheels for portability



Quick lock and release for handrails (optional)



Push wheelchair onto platform via two-way ramps

Features

- BMI
- Print
- Hold
- Tare
- Pre-Tare
- Backlight
- Portable Handle
- Wheel Castors
- Data Transfer

| Model | MS3830 Wheelchair Scale |
|---------------------|---|
| Capacity | 300 kg |
| Graduation | 100g |
| Overall Dimensions | 1150(D) x 800(W) x 66(H) mm |
| Platform Dimensions | 900(D) x 740(D) |
| Device Weight | 28.6 kg |
| Power Supply | AA batteries*6 / Power Adapter / Rechargeable battery pack (optional) |
| Accessories | Thermal Printer, Display Stand, Handrails |